

# Electrodermal Activity (EDA) Sensor Datasheet

EMG 03092020

## SPECIFICATIONS

- > **Range:** 0-25 $\mu$ S (@VCC=3V)
- > **Bandwidth:** 0-3Hz
- > **Consumption:**  $\pm$ 0.1mA
- > **Input Bias Current:**  $\pm$ 70pA
- > **CMRR:** 130dB
- > **Measurement:** continuous
- > **Current:** DC

## FEATURES

- > Skin resistance measurement
- > Unobtrusive & lightweight sensor
- > Pre-conditioned analog output
- > High signal-to-noise ratio
- > Ready-to-use form-factor
- > Medical-grade raw data output

## APPLICATIONS

- > Life sciences studies
- > Sympathetic nervous system monitoring
- > Human-Computer Interaction
- > Affective computing
- > Psychophysiology
- > Biomedical device prototyping
- > Arousal detection
- > Emotional cartography
- > Physiology studies
- > Relaxation biofeedback

## GENERAL DESCRIPTION

The biosignalsplux EDA sensor is capable of accurately measuring the electrical properties of the skin which changes. These changes are caused by alterations in sweat secretion and sweat gland activity as a result of changing sympathetic nervous system activity. The low-noise signal conditioning and amplification circuit design provide optimal performance in the detection of even the most feeble electrodermal skin response events.

## APPLICATION NOTES

The biosignalsplux EDA sensor is designed to acquire the change of skin activity such as sweat with two measuring electrodes. One example is the placement of the electrodes on the anterior side of the hand on two adjacent




Fig. 1. BIOSIGNALSPLUX EDA SENSOR (STANDARD VERSION).




Fig. 2. Typical raw EDA data (acquired with biosignalsplux).




Fig. 3. Example EDA placement on the index and middle finger.

**biosignalsplux**  
wearable body sensing platForm

PLUX – Wireless Biosignals, S.A.  
Av. 5 de Outubro, n. 70 – 2.  
1050-059 Lisbon, Portugal  
plux@plux.info  
<http://biosignalsplux.com/>

REV B

© 2020 PLUX

This information is provided "as is," and we make no express or implied warranties whatsoever with respect to functionality, operability, use, fitness for a particular purpose, or infringement of rights. We expressly disclaim any liability whatsoever for any direct, indirect, consequential, incidental or special damages, including, without limitation, lost revenues, lost profits, losses resulting from business interruption or loss of data, regardless of the form of action or legal theory under which the liability may be asserted, even if advised of the possibility of such damages.

# Electrodermal Activity (EDA) Sensor Data Sheet

EMG 0309202

fingers of interest (see Fig. 3). The electrodes are then connected to the sensor cable of the EDA sensor.

## TRANSFER FUNCTION

[0μS, 25μS]

$$EDA(\mu S) = \frac{\frac{ADC}{2^n} \times VCC}{0.12}$$

$$EDA(S) = EDA(\mu S) \times 1 \times 10^{-6}$$

VCC = 3V (operating voltage)

EDA(μS) – EDA value in microsiemens (μS)

EDA(S) – EDA value in in siemens (S)

ADC – Value sampled from the channel

n – Number of bits of the channel<sup>1</sup>

## ELECTRODE CONNECTIONS

| Electrode Cable | +   | -     |
|-----------------|-----|-------|
| Sleeve Color    | Red | Black |

## PHYSICAL CHARACTERISTICS

|                            |                                                         |                 |                 |
|----------------------------|---------------------------------------------------------|-----------------|-----------------|
| > W1 x L1 x H1:            | 1.6x2.2x0.5cm                                           | > W2 x L2 x H2: | 1.5x1.75x0.55cm |
| > W3:                      | 0.9cm                                                   | > L3:           | 0.5cm           |
| > A1:                      | 105.0±0.5cm                                             | > A2:           | 5.0±0.5cm       |
| > A3:                      | 5.0±0.5cm                                               | > D:            | 0.4cm           |
| > Available sleeve colors: | White, Black, Blue, Green, Red, Yellow, Gray, and Brown |                 |                 |
| > Weight:                  | 13g                                                     |                 |                 |

<sup>1</sup> The number of bits for each channel depends on the resolution of the Analog-to-Digital Converter (ADC); in biosignalsplux the default is 16-bit resolution (n=16), although 12-bit (n=12) and 8-bit (n=8) may also be found.